

User Charges Latest Addition:

| Sl. No | Instrument and Software | Total per sample (Rs) | |
|--------|--|-----------------------|---------|
| | | Academic | Private |
| 1 | Scanning Electron Microscope | 2000 | 3900 |
| 2 | Ultra-centrifuge | 2900 | 5800 |
| 3 | RNA sequencing | 5900 | 8260 |
| 4 | Exome sequencing | 5900 | 8260 |
| 5 | ChIP sequencing | 5900 | 8260 |
| 6 | sc-ATAC sequencing | 49000 | 51350 |
| 7 | sc-RNA sequencing | 49000 | 51350 |
| 8 | MD simulation | 9440 | 11800 |
| 9 | Server | 9440 | 11800 |
| 10 | GC-MS analysis | 5000 | 7000 |
| 11 | LC-MS/MS analysis, (Q-TRAP) | 5000 | 7000 |
| 12 | LC-MS/MS proteome analysis (digested sample) [ORBITRAP MS/MS] | 9000 | 12000 |
| 13 | LC-MS/MS proteome analysis (including sample preparation) [ORBITRAP MS/MS] | 12000 | 15000 |
| 14 | LC-MS/MS Quantitative proteomics (Labelled samples6 fractions maximum) [ORBITRAP MS/MS] | 35000 | 50000 |
| 15 | LC-MS/MS PTM analysis (Enriched sample) [ORBITRAP MS/MS] | 15000 | 20000 |
| 16 | LC-MS/MS Label free Quantitative proteome analysis (2sets of digested sample) [ORBITRAP MS/MS] | 20000 | 30000 |
| 17 | LC-MS/MS Label free Quantitative proteome analysis (including 2sets of sample preparation) [ORBITRAP MS/MS] | 25000 | 35000 |
| 18 | LC-MS/MS Label free Quantitative proteome analysis (3 sets of digested sample) [ORBITRAP MS/MS] | 30000 | 40000 |
| 19 | LC-MS/MS Label free Quantitative proteome analysis (including 3 sets of sample preparation) [ORBITRAP MS/MS] | 35000 | 45000 |
| 20 | LC-MS/MS Label free Quantitative proteome analysis (4 sets of digested sample) [ORBITRAP MS/MS] | 40000 | 50000 |
| 21 | LC-MS/MS Label free Quantitative proteome analysis (including 4 sets of sample preparation) [ORBITRAP MS/MS] | 45000 | 55000 |